



Accutome Connect™

DICOM Setup Quick Reference Guide

www.accutome.com

3222 Phoenixville Pike, Bldg 50 | Malvern, PA 19355 USA | info@accutome.com
800.979.2020 | 610.889.0200 | F: 610.889.3233 | European Office +31 (0) 485-350300



Setup DICOM Worklist and PACS

1. Launch Accutome Connect software.
2. On the Connect Home Screen select the Settings Button.



3. Press the “Facility Tab” and enter facility information.
(Name, Address, and Telephone Number)

A screenshot of the 'Facilities' tab in the software. It shows a form for entering facility information. The fields are: Name (ACCUTOME), Address (3222 Phoenixville Pike), City (Malvern), Zip (19355), State (PA), Country (USA), and Telephone (610 8890200).

4. Press “General” tab.
5. Select “DICOM Servers” tab.
6. The facility information that has been entered should show up on top.

A screenshot of the 'General' tab in the software. It shows the facility information entered in the previous step, including Name (ACCUTOME INC), Address (3222 Phoenixville Pike, Malvern, PA 19355, US), and Phone (800-979-2020). The 'PACS' tab is also visible.

7. Enter the PACS server information:
 - Local AE Title
 - Display as
 - IP Address or Host Name
 - AE Title
 - Port Number

(Information should be provided by your DICOM Administrator)

A screenshot of the 'PACS' tab in the software. It shows a form for entering PACS server information. The fields are: Local AE Title (ConnectUnit), Remote Servers (List of Remote Servers), Display as (Connect), IP address or host name (10.108.1.14), AE title (ConnectUnit), and Port number (61200). There are buttons for 'New Server', 'Save Server', 'Remove Sel', 'Set Selected as Default', and 'Verify Connection'. The DICOM Folder is set to C:\Accutome\DICOM\.

8. If both PACS and Worklist are using the same server, you will be able to copy the server information by selecting “Duplicate to Worklist Server”.



9. Select "Save Server"


- PACS server information will copy into the Worklist server information fields.
- Press "Set Selected as Default" button.
- Press "Verify Connection" button.

(If "Duplicate to Worklist Server" was selected, information will be copied in worklist fields already. Make sure you change the "Port Number" to the correct worklist port number)



You should see a pop-up window that says Success


If you see any other message, check that your Server information is entered correctly.



10. If PACS and Worklist are on two separate servers,

- Enter Worklist Server Information:
 - Local AE Title
 - Display As
 - IP Address or Host Name
 - AE Title
 - Port Number


(Information should be provided by your DICOM Administrator)



11. Press "Save" button.


12. Press "Set Select as Default" button.

13. Press "Verify Connection" button.



You should see a pop-up window that says Success

If you see any other message, check that your Server information is entered correctly.



14. Press "Done" to save changes and exit setup.



For assistance with the device upgrade or any other questions please feel free to contact Accutome diagnostic support at diagnostics@accutome.com or 610-889-0200.