

UBM Recommended Computer Specs

1. **Processor**

Intel Core 2 Duo 2.20 GHz, Core i3 2.2 GHz or higher

2. **Memory**

3.0 GB, DDR2-667 SDRAM or higher

3. **Video Card or Chip**

Desktop PC: Any NVIDIA or ATI Video Card with 256MB graphic memory or more.

Laptop PC: Any NVIDIA, ATI or Intel Video Chip

4. **Screen Resolution of Monitor**

High resolution monitor (minimum 1366 X 768 pixels; recommended 1600 X 900 pixels or higher) 19 inch for desktop
15.4 inches or higher screen size for laptop.

5. **Hard Disk**

Minimum 80 GB, 7200 RPM, recommended 160 GB or higher for high patient flow site (more than 30 B-Scans per month).

6. **USB Ports**

Minimum of 2 USB 2.0

7. **Medical Grade Power Supply**

8. **Software and Operation System (OS)**

Microsoft Word[®] 2003, 2007 and 2010 (for report generation)

Windows 7[®], Microsoft Windows[®] Vista or Microsoft Windows[®] XP Professional Edition, SP2, Microsoft Windows[®] 8.

UBM Recommended Computer Specs

If the computer on which you intend to load the UBM Software does not meet or exceed the listed specifications, we highly recommend that you do not install the UBM plus Software.

It is recommended that the computer used in conjunction with the UBM Plus system be dedicated for the sole use of Accutome's UBM. Installation of software not vital to the operation of the UBM is highly discouraged, as it may detrimentally affect the performance of the UBM Plus system. Additionally, accessing the Internet will greatly affect the performance of the UBM and open the system to viruses and other harmful external factors and is therefore not recommended by Accutome.