

# ClearScan® Probe Cover

Item # 24-6350

ClearScan® is an innovative single use ultrasound probe cover which eliminates the difficulties and problems of the shell and gel system. The ClearScan® consists of an extremely thin film that is acoustically invisible, resulting in the same distortion free ultrasound image as with the shell/gel technique, but without the discomfort. ClearScan® also has the ability to scan 10mm posterior from the limbus.

## Open Shell Challenges:

- Patient must be in reclined position
- Cannot be used on the eye directly after surgery
- Requires BSS/distilled water to be used in shell
- Shell must be inserted onto patient eye causing discomfort

## ClearScan Solutions:

- Sterile probe cover can help to reduce the risk of patient to patient transfer of infection
- Internal bag pressure can be easily reduced for eyes with low IOP
- Examine delicate structures following recent surgery
- Examine young children
- Simple exam setup, approximately 30 seconds

Distributed by:

