



Get a Grip on Trephines!

Accutome would like to announce that we now carry Surgistar's line of radial vacuum trephines and punches. These vacuum trephines consistently create accurate and uniform cuts by dividing the corneal chamber into 16 chambers, thus reducing corneal distortion.

Features:

- Ergonomic grip for easier control in precision cuts on patient corneas
- 16 radial marks for easy suture placement
- Fine cross hairs for easy, accurate alignment on visual axis
- Increased blade size font for easy size recognition

Single-use, sterile trephines and punches are available in an assortment of sizes.

Visit www.accutome.com to order online or call your Accutome representative now.



1-800-979-2020 • 610-889-0200
+31 [0] 485-350300 Europe
FAX 610-889-3233 • info@accutome.com