


AccuDrop

Introducing Accutome's Dropless Therapy Cannula

Stop, **Drop** & *Go!*



Visit  **ACCUTOME**[®]
A HALMA COMPANY
At Booth Number 3144

AccuDrop

Stop, Drop & Go!

- 25 gauge* blunt tip cannula
- Side port opening allows for gentle drug dispersion
- Blunt tip design to reduce zonular trauma
- Angled shaft designed for targeted drug delivery
- 90° bend facilitates easier insertion into anterior chamber
- Designed by Accutome for administration of dropless cataract drug formulation

**Other sizes available*



ACCUTOME[®]
A HALMA COMPANY

www.accutome.com

800.979.2020 • 610.889.0200

F:610.889.3233 • info@accutome.com

Request A Sample

Name: _____

Company: _____

E-mail: _____

Phone: _____

Address: _____



LIKE US

facebook.com/pages/Accutome-Inc



FOLLOW US

linkedin.com/company/accutome-inc



WATCH RESOURCES

youtube.com/user/Accutome



TWEET US

twitter.com/accutome