

AccuPen[®]

Handheld Tonometer

Item # 24-3000

Features:

- Gravity Offset Technology provides accurate IOP measurements with less calibration in any testing position.
- New ergonomic design allows for easy visualization of the cornea.
- Low cost battery takes up to 15,000 readings.
- Self Calibrating. No more up down and around.
- Large LCD Display stores up to 9 measurements with a calculated average.
- Adjusted IOP Calculation to correct for varying CCTs.
- Comes with heavy duty traveling case and lanyard for easy transportation.

Specifications:

- Weight: 3.0oz.
- Physical Size: 7 ¼" x 1 ¼"
- Accuracy: +/- 2mmHg [compared to Goldman Tonometer]*

*When using an average of readings.

The AccuPen is supported by Accutome's Superior Customer Service and Technical Experience. Call to order the AccuPen today and receive 100 free AccuTip covers.



www.accutome.com

1-800-979-2020 • 610-889-0200 • +31 (0) 485-350300 • FAX 610-889-3233 • info@accutome.com

ACCUTOME[®]
A HALMA COMPANY