



The Complete Biometry Solution

Seamless workflow integration with the IOLMaster 500 from ZEISS

A-Scan Plus Connect®

Item Number: 24-4400

Help to ensure accuracy¹ and improve efficiency with the new A-Scan Plus Connect®. The Connect ultrasound biometer provides a data link to the ZEISS IOLMaster® 500 for seamless workflow integration. The universal ultrasound interface of the IOLMaster allows for quick and easy data transfer between the IOLMaster 500 and the A-Scan Plus Connect.

Features:

- Immersion and contact modes.
- 100% Measurement capability regardless of lens opacities or fixation difficulties.
- Industry leading resolution.
- User-friendly interface.
- Less patient chair time – faster measurement capture².
- Automatic alignment detection and sclera recognition eliminate marginally aligned scans.
- Portable, lightweight design – plug into any Windows based pc, laptop or tablet device.
- Share information easily – adaptable document transfer via EMR, email or printer.
- Modern third and fourth generation formulas including the Hoffer®Q, SRK/T, Holladay, Haigis as well as new post refractive formulas.
- Optimized lens constants for superior surgical outcomes.
- Unlimited software licenses.
- Streamline your patient data - shares information with IOLMaster 500, B-Scan Plus and UBM Plus.

¹Haigis W, Mlynski J, Comparative axial length measurements using optical and acoustic biometry in normal persons and in patients with retinal lesions, White Paper, Carl Zeiss Meditec, 2009

²When compared to original A-Scan Plus

3222 Phoenixville Pike • Malvern PA 19355 USA
800.979.2020 • 610.889.0200 • F:610.889.3233 • www.accutome.com

 **ACCUTOME**[®]
A HALMA COMPANY

EN_32_025_0090I