



## The Accutome Connect Platform is DICOM Compatible.

The DICOM (Digital Imaging and Communications in Medicine) is the standard for communication of medical imaging devices, by specifying a protocol for network communications, media communications, file formats and other information. Accutome's A-Scan Plus Connect, UBM Plus, and B-Scan Plus has DICOM compatible software.

### DICOM Capabilities

- Export and store to your DICOM server or PACS/ image management system.
- Verification of connection to other Application Entities (AEs).
- The only A-Scan to export axial length and reports to the DICOM/PACS, allowing manipulation of data and calculations within your DICOM system.

### Features

- DICOM exporting capacity
- One license one location. No need for multiple license purchases
- Easy installation. No need to send unit in for upgrade. Can be completed remotely.
- Easy activation. Minimal information needed to activate exporting capabilities.
- 100% product support via Accutome home office.

***As always, our product specialists are here to support your installation and training.***

3222 Phoenixville Pike • Malvern PA 19355 USA • 800.979.2020  
610.889.0200 • F:610.889.3233 • [www.accutome.com](http://www.accutome.com) • [info@accutome.com](mailto:info@accutome.com)

 **ACCUTOME**<sup>®</sup>  
A HALMA COMPANY