

## **B-Scan Plus**

### **Recommended Computer Specs**

1. **Processor**

Intel Core 2 Duo 2.20 GHz, Core i3 2.2 GHz or higher

2. **Memory**

3.0 GB, DDR2-667 SDRAM or higher

3. **Video Card or Chip**

Desktop PC: Any NVIDIA or ATI Video Card with 256MB graphic memory or more.

Laptop PC: Any NVIDIA, ATI or Intel Video Chip

4. **Screen Resolution of Monitor**

High resolution monitor (minimum 1366 X 768 pixels; recommended 1600 X 900 pixels or higher) 19 inch for desktop  
15.4 inches or higher screen size for laptop.

5. **Hard Disk**

Minimum 80 GB, 7200 RPM, recommended 160 GB or higher for high patient flow site (more than 30 B-Scans per month).

6. **USB Ports**

Minimum of 2 USB 2.0

7. **Medical Grade Power Supply**

8. **Software and Operation System (OS)**

Microsoft Word<sup>®</sup> 2003, 2007 and 2010 (for report generation)  
Windows 7<sup>®</sup>, Microsoft Windows<sup>®</sup> Vista or Microsoft Windows<sup>®</sup> XP Professional Edition, SP2, Microsoft Windows<sup>®</sup> 8.

## **B-Scan Plus Recommended Computer Specs**

**If the computer on which you intend to load the B-Scan Software does not meet or exceed the listed specifications, we highly recommend that you do not install the B-Scan Plus Software.**

**It is recommended that the computer used in conjunction with the B-Scan Plus system be dedicated for the sole use of Accutome's B-Scan. Installation of software not vital to the operation of the B-Scan is highly discouraged, as it may detrimentally affect the performance of the B-Scan Plus system. Additionally, accessing the Internet will greatly affect the performance of the B-Scan and open the system to viruses and other harmful external factors and is therefore not recommended by Accutome.**