# B-Scan Vet Recommended Computer Specs

#### 1. Processor:

Intel<sup>®</sup>Core<sup>TM</sup> 2 Duo 2.00 GHz or higher

#### 2. Memory:

2.0 GB, DDR2 SDRAM

#### 3. Video:

256 MB NVIDIA GeForce 8600M GT graphics card WSXGA (1280 X 800 pixels minimum) monitor 15.4 inches or higher screen size.

### 4. Hard Disk:

Minimum 80 GB, 7200 RPM, recommended 160 GB or higher for high patient flow site (more than 30 B-Scans per month).

## 5. Optical Drive:

8X (or higher) DVD+/-RW with control software (needed for archiving)

### 6. **USB Ports**

Minimum of 3 USB 2.0

# 7. Medical Grade Power Supply

#### 8. Software

Microsoft Word<sup>®</sup> 2003 or 2007 (needed for report generation) Microsoft Windows<sup>®</sup> XP Professional Edition, SP2 Microsoft Windows<sup>®</sup> Vista Business or Ultimate editions Microsoft Windows 7<sup>®</sup>

### 9. File System

NTFS file System

\*Does NOT support either Dual Monitors or Dual Video Cards

# B-Scan Vet Recommended Computer Specs

If the computer on which you intend to load the B-Scan Software does not meet or exceed the listed specifications, we highly recommend that you do not install the B-Scan Vet Software.

If you wish to purchase an Accutome recommended Dell Latitude E6500 laptop, please contact Dell at 1-800-915-3355

Accutome recommends using the Xerox 8560 and Lexmark D830 printer.

If you wish to purchase an Accutome recommended Printer 8560, please contact Drew Armer at CDW at (877) 685-8896 and reference the "Accutome Xerox Printer."

The computer used in conjunction with the B-Scan Vet system should be dedicated for the sole use of Accutome's B-Scan. Installation of software not vital to the operation of the B-Scan is highly discouraged, as it may detrimentally affect the performance of the B-Scan Vet system. Additionally, accessing the Internet will greatly affect the performance of the B-Scan and open the system to viruses and other harmful external factors and is therefore not recommended by Accutome.