



## PSL Classic® Portable Slit Lamp

Item Number: 24-1450

### Features:

- *Precision machined aluminum chassis*  
Designed and manufactured with the strength to withstand the daily traumas of the demanding ophthalmic environment. Lightweight, durable and robust.
- *Advanced optics, x10 & x16 magnifications*  
Top quality multi layered glass is crafted into our multi-element and multi-layer coated optical system to give maximum performance. Clarity is guaranteed.
- *Controllable illumination*  
Proven halogen illumination system with over 2,000 hours of bulb life; illumination levels are controlled by a rheostat giving you control of your examination parameters. From zero to maximum, you choose.
- *Big Slit Lamp features, portable usability*  
Unique features such as the fixation targets, and the 1mm square light patch for assessing a/c flare set the new standard for Portable Slit Lamps.