

Tonoclear

STERILE




Part Number: 2414-P-7045

Product Overview:


A single-use tonometer prism to be used in combination with a GAT, to determine the IOP. The Tonoclear prism is designed to accommodate variations in light intensity from the slit lamp to ensure a clear view of the mires at all times. The individually packed prisms come pre-assembled to reduce any risk of a cross-contamination and hassle in your workflow.

Features:

- Single-use prisms (box of 100)
- Pre-assembled
- Clear visualization of the mires
- Attenuating sleeve to reduce sensitivity to light variations
- Individually packaged and sterilized with Ethylene Oxide
- Sterilization validity for 3 years
- Compatible with all GATs



Tonoclear mires image



Competitor mires

For more information, visit www.keelerusa.com, call us at 800-523-5620, or email us at customerservice@keelerusa.com