

PachPen[®]

Handheld Pachymeter

Item # 24-5100

Accutome's PachPen is the perfect combination of accuracy, portability and comfort. The superior ergonomic design fits comfortably in any hand, which is a major improvement from other handheld models. When you consider the PachPen is equipped with Accutome's Digital Signal Analysis, there is no reason to use anything else at your facility.

Features:

- Accurate - 65 MHz sampling probe ensures precise measurements.
- IOP Correction Calculation - Converts IOP measurements in seconds.
- Gentle Touch - 2.5mm diameter probe helps to ensure patients comfort.
- No need to recharge - Long lasting lithium battery lasts approximately 15,000 readings.
- Lightweight - The PachPen weighs only 3.0 oz, the lightest on the market.

